

Oxygen Mask with Reservoir Bag

PRODUCT DATA SHEET



DESCRIPTION

m|devices non-breathing Oxygen Mask with Reservoir Bag is available with 2.1 m kink resistant star lumen and is made from non-toxic medical grade PVC. The mask is elongated in shape and provides higher concentration of oxygen than is able to be provided with a standard face mask or nasal prongs and is available in two sizes; paediatric and adult.

FEATURES AND BENEFITS

- Paediatric and adult version available
- Adjustable straps
- Removable vent covers x 2
- Inlet valve: non-return
- With reservoir bag

PRODUCT DATA

- Tubing length: 2.1m
- Tubing diameter: 3.94cm
- Tubing type: Star lumen tubing
- Material: PVC
- DEHP: Yes
- Latex: No
- Shelf Life: 5 Years
- ARTG Number: 228231

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	TUBING LENGTH	SIZE	CARTON/QTY
AN62001NS	Oxygen Mask with Reservoir Bag Elongated Shape with Tubing	2.1m	Paediatric	50
AN62003NS	Oxygen Mask with Reservoir Bag Elongated Shape with Tubing	2.1m	Adult	50