


Nasal Oxygen Cannulas

PRODUCT DATA SHEET



DESCRIPTION

m|devices Nasal Oxygen Cannulas are made from non-toxic medical grade PVC and are designed to simplify oxygen administration and are super soft and flexible to maximise patient comfort. The green coloured, kink resistant star lumen tubing increases visibility and can assist in reducing connection errors.



FEATURES AND BENEFITS

- 3 sizes available; Infant, Child and Adult
- Super soft, flexible nasal prong inserts for added patient comfort
- Kink resistant star lumen tubing

PRODUCT DATA

- Tubing length: 2.1m, 3m
- Tubing type: Star Lumen
- Distance between nasal prongs and internal diameter (ID) of nasal prong tip:
INFANT 8mm, ID 1.6mm
CHILD 12mm, ID 2.0mm
ADULT 16mm, ID 3.0mm
- DEHP: Yes
- Latex: No
- Shelf Life: 5 Years
- ARTG Number: 228231

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	TUBING LENGTH	SIZE	CARTON/QT
AN070000NS	Nasal Oxygen Cannula	2.1m	Infant	50
AN070001NS	Nasal Oxygen Cannula	2.1m	Child	50
AN070002NS	Nasal Oxygen Cannula	2.1m	Adult	50
AN070003NS	Nasal Oxygen Cannula	3m	Infant	50
AN070004NS	Nasal Oxygen Cannula	3m	Child	50
AN070005NS	Nasal Oxygen Cannula	3m	Adult	50