

Nasal Oxygen Cannulas

DESCRIPTION

m|devices Nasal Oxygen Cannulas are made from non-toxic medical grade PVC and are designed to simplify oxygen administration and are super soft and flexible to maximise patient comfort.

FEATURES & BENEFITS

- 3 sizes available; Infant, Child and Adult
- Super soft, flexible nasal prong inserts for added patient comfort
- Kink resistant star lumen tubing


PRODUCT DATA

- Tubing length: 2.1m, 3m
- Tubing type: Star Lumen
- Distance between nasal prongs and internal diameter (ID) of nasal prong tip:
INFANT 8mm, ID 1.6mm
CHILD 12mm, ID 2.0mm
ADULT 16mm, ID 3.0mm
- DEHP YES = Product contains Phthalates such as the PVC Plasticizer Di-(2-ethylhexyl)-phthalate (DEHP).
- Not Made with Natural Rubber Latex
- ARTG Number: 228233
- Shelf Life: 5 Years

STANDARDS

- EN ISO 11135-1:2007 Sterilization of Health Care Product
- EN ISO11607-1:2006 Packaging for terminally sterilized medical devices — Part 1: Requirements for materials, sterile barrier systems and packaging systems
- ISO 10993-1:2009 - Biological evaluation of medical devices

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Nasal Tip Distance 8mm-16mm
2	Size Adjuster
3	Connector
4	Tube

Nasal Oxygen Cannulas

Nasal Oxygen Cannulas				
RE-ORDER	DESCRIPTION	TUBING LENGTH	SIZE	CARTON/QTY
AN070000NS	Nasal Oxygen Cannula	2.1m	Infant	50
AN070001NS	Nasal Oxygen Cannula	2.1m	Child	50
AN070002NS	Nasal Oxygen Cannula	2.1m	Adult	50
AN070003NS	Nasal Oxygen Cannula	3m	Infant	50
AN070004NS	Nasal Oxygen Cannula	3m	Child	50
AN070005NS	Nasal Oxygen Cannula	3m	Adult	50