

# Nebuliser Kits

## DESCRIPTION

m|devices Nebuliser Kits are available in 4 different mask sizes; Paediatric Standard Shape, Paediatric Elongated Shape, Adolescent Standard shape and Adult Elongated Shape. A mouth piece only kit is also available. They all include 2.1m and 3m kink resistant star lumen tubing and a 10mL jar. All kits are made from medical grade PVC.

## FEATURES & BENEFITS

- 4 sizes available ranging from Paediatric to Adult
- Standard shape and Elongated shape options
- 10mL jar included
- Kink resistant star lumen tubing
- Mouth piece versions available AN082000NS
- Standard connectors


## PRODUCT DATA

- Tubing length: 2.1m, 3m
- Tubing type: Star Lumen
- Tubing diameter: 3.94mm
- DEHP YES = Product contains Phthalates such as the PVC Plasticizer Di-(2-ethylhexyl)-phthalate (DEHP).
- Not Made with Natural Rubber Latex
- ARTG Number: 228231
- Shelf Life: 5 Years
- Material: PVC

## STANDARDS


- EN 980:2008 Symbols for use in the labelling of medical devices
- ISO 10993-1:2003 - Biological evaluation of medical devices
- ISO 10993-5:1999 - Biological evaluation of medical devices Part 5: Tests for in vitro cytotoxicity.

## PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Mask
2	Nose Clip
3	Elastic Strap
4	Nebuliser Jar
5	Tubing
6	Connector
7	Connector for Mask

Nebuliser Kit with Mask



ITEM	DESCRIPTION
1	Mouthpiece
2	Nebuliser Jar
3	Tubing
4	Connector

Nebuliser Kit with Mouthpiece

# Nebuliser Kits

Nebuliser Kits				
RE-ORDER	DESCRIPTION	TUBING LENGTH	SIZE	CARTON/QTY
<b>AN080000NS</b>	Nebuliser Kit - Standard Shape 10mL NebuliserJar	2.1m	Paediatric	50
<b>AN080001NS</b>	Nebuliser Kit - Elongated Shape 10mL NebuliserJar	2.1m	Paediatric	50
<b>AN081000NS</b>	Nebuliser Kit - Standard Shape 10mL NebuliserJar	2.1m	Adolescent	50
<b>AN081001NS</b>	Nebuliser Kit - Elongated Shape 10mL NebuliserJar	2.1m	Adult	50
<b>AN082000NS</b>	Nebuliser Kit - Mouth Piece 10mL NebuliserJar	2.1m	-	50
<b>AN083000NS</b>	Nebuliser Kit - Standard Shape 10mL NebuliserJar	3m	Paediatric	50
<b>AN083001NS</b>	Nebuliser Kit - Elongated Shape 10mL NebuliserJar	3m	Paediatric	50
<b>AN084000NS</b>	Nebuliser Kit - Standard Shape 10mL NebuliserJar	3m	Adolescent	50
<b>AN084001NS</b>	Nebuliser Kit - Elongated Shape 10mL NebuliserJar	3m	Adult	50
<b>AN085000NS</b>	Nebuliser Kit - Mouth Piece 10mL NebuliserJar	3m	-	50