

Guedel Airways

PRODUCT DATA SHEET



DESCRIPTION

The m|devices Guedel Airways are made from semi-rigid, non-toxic polyethylene and are compliant to ISO standard lengths. They are available in a wide range of sizes which are colour coded for easy identification, and are suitable for both adult and paediatric patients.



FEATURES AND BENEFITS

- Suitable for adult and paediatric patients
- Semi-rigid, non toxic polyethylene
- Colour coded
- Wide range of sizes

PRODUCT DATA

- DEHP: Yes
- Latex free
- Sterile: No
- Material: Polyethylene
- Shelf Life: 5 Years
- ARTG Number: 227534

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	COLOUR	LENGTH	CARTON/QTY
AN090000NS	Guedel Airway	No.000	Pink	40mm	50
AN090001NS	Guedel Airway	No.00	Blue	50mm	50
AN090002NS	Guedel Airway	No.0	Black	60mm	50
AN090003NS	Guedel Airway	No.1	White	70mm	50
AN090004NS	Guedel Airway	No.2	Green	80mm	50
AN090005NS	Guedel Airway	No.3	Yellow	90mm	50
AN090006NS	Guedel Airway	No.4	Red	110mm	50
AN090007NS	Guedel Airway	No.5	Light Blue	120mm	50
AN090008NS	Guedel Airway	No.6	Orange	130mm	50