

# Intubating Stylets

## PRODUCT DATA SHEET



### DESCRIPTION

m|devices Intubating Stylets are available in three different sizes, 6Fr, 10Fr and 14Fr to accommodate a wide range of endotracheal tube sizes. They are made of malleable aluminium covered in PVC sheath and are easily shaped to assist with difficult intubations.


### FEATURES AND BENEFITS

- Available in 3 sizes: 6Fr, 10Fr, 14Fr
- Made of malleable aluminium with PVC sheath
- Easy to bend into desired curve required to assist difficult intubation
- Suitable for a wide range of ET tubes:
  - AN100000 6Fr for ETT 3.0mm - 4.0mm
  - AN100001 10Fr for ETT 4.5mm - 7.5mm
  - AN100002 14Fr for ETT 8.0mm - 10.0mm

### PRODUCT DATA

- Lengths: 30cm or 38cm
- DEHP: Yes
- Latex free
- Sterile: EO
- Material: Aluminium, PVC
- Shelf Life: 5 Years
- ARTG: 336796

### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
AN100000	Intubating Stylet for ETT 3.0mm - 4.0mm	6Fr	20
AN100001	Intubating Stylet for ETT 4.5mm - 7.5mm	10Fr	20
AN100002	Intubating Stylet for ETT 8.0mm - 10.0mm	14Fr	20