

Intubating Stylets

DESCRIPTION

m|devices Intubating Stylets are available in three different sizes, 6Fr, 10Fr and 14Fr to accommodate a wide range of endotracheal tube sizes. They are made of malleable aluminium covered in PVC sheath and are easily shaped to assist with difficult intubations.

FEATURES & BENEFITS

- Available in 3 sizes: 6Fr, 10Fr, 14Fr
- Made of malleable aluminium with PVC sheath
- Easy to bend into desired curve required to assist difficult intubation
- Suitable for a wide range of ET tubes:
 - AN100000 6Fr for ETT 3.0mm - 4.0mm
 - AN100001 10Fr for ETT 4.5mm - 7.5mm
 - AN100002 14Fr for ETT 8.0mm - 10.0mm


PRODUCT DATA

- Lengths: 30cm or 38cm
- DEHP YES = Product contains Phthalates such as the PVC Plasticizer Di-(2-ethylhexyl)-phthalate (DEHP).
- Not Made with Natural Rubber Latex
- ARTG Number: 336796
- Shelf Life: 5 Years
- Material: Aluminium, PVC
- Sterile: EO

STANDARDS

- EN ISO 14155-1:2003 Clinical investigation for medical devices for human subjects – part1 General requirements
- ISO 14971:2007 Medical devices – Application of risk management to medical devices - Part 1: Evaluation and testing
- ISO 10993-1:2003 Biological evaluation of the medical devices – Part 1: Guidance on selection of tests

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Malleable Tip

Intubating Stylets

Intubating Stylets			
RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
AN100000	Intubating Stylet for ETT 3.0mm - 4.0mm	6Fr	20
AN100001	Intubating Stylet for ETT 4.5mm - 7.5mm	10Fr	20
AN100002	Intubating Stylet for ETT 8.0mm - 10.0mm	14Fr	20