

Oxygen Tubing



PRODUCT DATA SHEET

DESCRIPTION

M|devices oxygen tubing is available in three different lengths and are made from a green non-toxic medical grade PVC, is kink resistant star lumen tubing. The green colour increases visibility and assists in reducing connection errors. The 3m length also comes with an optional step connector.


FEATURES AND BENEFITS

- Three lengths available
- Option step connector for 3 m length
- Kink resistant star lumen tubing
- Soft tubing
- Medical grade PVC

PRODUCT DATA

- Tubing length: 2.1, 3 or 4m
- Material: PVC
- DEHP: Yes
- Shelf life: 5 years
- ARTG Number: 267877

PRODUCT PACKAGING LABEL



RE-ORDER	SIZE	CARTON QTY
AN130000NS	2.1m	50pcs
AN130001NS	3m	50pcs
AN130004NS	4m	50pcs