

Oxygen Tubing

Non Sterile

DESCRIPTION

m|devices Oxygen Tubing is available in three different lengths, are made from green, non-toxic, medical grade PVC, and is kink resistant star lumen tubing. The 3m length also comes with an optional step connector.

FEATURES & BENEFITS

- Three lengths available
- Optional step connector for 3m length
- Soft tubing
- Kink resistant star lumen tubing
- Medical grade PVC


PRODUCT DATA

- Tubing length: 2.1m, 3m, 4m, 30m
- Tubing type: Star Lumen
- Material: PVC
- ARTG Number: 267877
- Shelf Life: 5 Years
- DEHP YES = Product contains Phthalates such as the PVC Plasticizer Di-(2-ethylhexyl)-phthalate (DEHP).

STANDARDS

- EN 980:2008 Symbols for use in the labelling of medical devices
- EN ISO 14644-1:1999 Cleanrooms and associated controlled environments — Part 1: Classification of air cleanliness.
- ISO 14971:2007 Medical devices — Application of risk management to medical devices

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Tubing
2	Connector

Oxygen Tubing

Non Sterile

Oxygen Tubing			
RE-ORDER	DESCRIPTION	TUBING LENGTH	CARTON/QTY
AN130000NS	Oxygen Tubing	2.1m	50
AN130001NS	Oxygen Tubing	3m	50
AN130004NS	Oxygen Tubing	4m	50
AN130030NS	Oxygen Tubing	30m	4
AN135001NS	Oxygen Tubing with White Step Connector	3m	50