

Bubble Oxygen Tubing Green

DESCRIPTION

m|devices Bubble Oxygen Tubing has incremental segmentations of two diameters (3.5mm and 7.5mm) throughout the length of tubing at 1m increments. This is specifically designed to enable the clinician to cut the tubing at the desired length and tubing diameter required to fit with the oxygen mask/nebuliser of choice and the corresponding oxygen source.

FEATURES & BENEFITS

- Flexible tubing that can be easily cut at desired length
- Allows complete customisation to connect to selected product of choice
- Coiled in dispenser box to assist with ease of removal


PRODUCT DATA

- Oxygen tubing length: 30m
- Segment diameter: ID 3.5mm - 7mm OD 5.3mm - 1.5mm
- Non Sterile
- Material: PVC
- Not made with Natural Rubber Latex
- Shelf Life: 5 Years
- ARTG Number: 267877

STANDARDS

- ISO 14971:2007 Medical devices-Application of risk management to medical devices
- ISO10993-1:2003 Biological evaluation of the medical devices-Part 1: Guidance on selection of tests
- ISO10993-10:2013 Biological evaluation of medical devices Part 10: Tests for irritation and skin sensitisation

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	ID - Small 3.5mm
2	ID - Large 7.0mm

Bubble Oxygen Tubing

RE-ORDER	DESCRIPTION	LENGTH	COLOUR	CARTON/QTY
AN131030NS	Bubble Oxygen Tubing	30m	Green	4