

Bubble Oxygen Tubing

PRODUCT DATA SHEET



DESCRIPTION

m|devices Bubble Oxygen Tubing has incremental segmentations of two diameters (3.5mm and 7mm) throughout the length of the tubing at 1m increments. This is specifically designed to enable the clinician to cut the tubing at the desired length and tubing diameter required to fit with the oxygen mask/nebulizer of choice and the corresponding oxygen source.


FEATURES AND BENEFITS

- Flexible tubing that can be easily cut at desired length
- Allows complete customisation to connect to selected product of choice
- Coiled in dispenser box to assist with ease of removal
- Latex free

PRODUCT DATA

- Oxygen tubing length: 30m
- Segment diameters: ID 3.5mm – 7mm
OD 5.3mm- 0.5mm
- Material: PVC
- DEHP: Yes
- Shelf Life: 3 years
- ARTG number: 228233

PRODUCT PACKAGING LABEL



RE-ORDER	SIZE	CARTON QTY
AN131030NS	30m	4 pcs