

Low Profile Anaesthetic Masks

PRODUCT DATA SHEET



DESCRIPTION

m|devices Low Profile Anaesthetic Masks are colour coded for easy identification. They are transparent, slim line and moulds to the face. This unique design offers an anaesthetic mask that is ergonomically shaped with a wide cuff that contours to the patient's face to obtain an effective seal. Unlike traditional cushioned masks, this range is not dependant upon the cushion to have the adequate air volume provided and eliminates the need to adjust the cushion by inflating or deflating prior or during use.

FEATURES AND BENEFITS

- Low profile mask
- Cushion-less design requires no inflation/deflation of cushion prior or during use
- Ergonomically shaped to contour to the patient's face to obtain an effective seal
- Available in 4 sizes to cover all face sizes and shapes
- Neutral scent
- Colour coded
- Easy grip design

PRODUCT DATA

- PVC free
- Latex free
- DEHP free
- Connector size:
Paediatric 15mm
Child, Small Adult and Large Adult 22mm
- Complies with ISO 5365-1
- Material: Medical Grade PP & TPE
- Shelf Life: 3 Years
- ARTG Number: 228655
- Non-Sterile

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	COLOUR	CARTON/QTY
AN145001NS	Low Profile Anaesthetic Mask	Paediatric	Pink	30
AN145002NS	Low Profile Anaesthetic Mask	Child	Green	30
AN145003NS	Low Profile Anaesthetic Mask	Small Adult	Purple	35
AN145004NS	Low Profile Anaesthetic Mask	Large Adult	Grey	35