

Anaesthetic Circuit Bain/Coaxial Circuit with T Connector

PRODUCT DATA SHEET



DESCRIPTION

The m|devices anaesthetics circuits are available in a range of different configurations to meet clinical practical needs. The collapsible bain/coaxial circuits come with a T-connector and a 2L bag and three differing tubing lengths to meet all your clinical needs.


FEATURES AND BENEFITS

- Three lengths available: 1.8m, 2.4m and 3m
- T connector to 60 cm collapsible circuit
- 2L bag
- Latex free

PRODUCT DATA

- | | | | |
|-------------|---------|----------------|---------|
| ■ Material: | PE, EVA | ■ Shelf life: | 5 Years |
| ■ DEHP: | Free | ■ ARTG Number: | 269092 |

PRODUCT PACKAGING LABEL



RE-ORDER | SIZE | CARTON QTY

AN150003NS	1.8m	20pcs
AN150004NS	2.4m	20pcs
AN150005NS	3m	20pcs
AN151005NS	3m bag limb	10pcs