

Anaesthetic Circuit

Bain/Coaxial Circuit with T Connector

DESCRIPTION

m|devices Anaesthetic Circuits are available in a range of different configurations to meet clinical practical needs. The collapsible Bain/Coaxial circuits with a T Connector and a large 2L bag with different tubing lengths to meet all your clinical needs.


FEATURES & BENEFITS

- Three lengths available: 1.8m, 2.4m and 3m
- T Connector to 60cm collapsible circuit
- 2L bag

PRODUCT DATA

- Not made with DEHP
- Shelf Life: 3 Years
- Material: PE, EVA
- ARTG Number: 269092

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	T-Connector
2	Bain Coaxial Circuit 1.8m
3	Limb
4	Reservoir Bag
5	Elbow connector with Port Cap Luer

Anaesthetic Circuit

RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
AN150003NS	Anaesthetic Circuit Bain/Coaxial Circuit with T Connector	1.8m	20
AN150004NS	Anaesthetic Circuit Bain/Coaxial Circuit with T Connector	2.4m	20
AN150005NS	Anaesthetic Circuit Bain/Coaxial Circuit with T Connector	3m	20
AN151005NS	3m Bag Limb - White	3m	10