

Anaesthetic Circuit

Corrugated Circuit with Rotatable Y Connector

PRODUCT DATA SHEET



DESCRIPTION

The m|devices anaesthetics circuits are available in a range of different configurations to meet clinical practical needs, offering adult and paediatric solutions. The collapsible corrugated circuits come with a rotating wye connector, available in three lengths and come with a 2L bag.


FEATURES AND BENEFITS

- Three lengths available for adult circuit: 1.8m, 2.4m and 3m
- Tube length for paediatric circuit: 1.8m
- Rotatable wye connector
- Latex free

PRODUCT DATA

- | | | | |
|------------------|----------------|----------------|---------|
| ▪ Tubing length: | 1.8, 2.4 or 3m | ▪ Shelf life: | 5 Years |
| ▪ Material: | PE, EVA | ▪ ARTG Number: | 269092 |
| ▪ DEHP: | Free | | |

PRODUCT PACKAGING LABEL



RE-ORDER | SIZE | CARTON QTY

AN150006NS	1.8m	20pcs
AN150007NS	2.4m	10pcs
AN150012NS	3m	10pcs
AN151005NS	1.8m	10pcs