

Yankauer Suction Handles Straight Tip

DESCRIPTION

m|devices Yankauer Suction Handles consists of suction tube. The suction tips are available as a straight tip offering with vent, without vent, or with an on/off switch. Smooth and clean edges for atraumatic suctioning made from a transparent medical grade PVC.

FEATURES & BENEFITS

- Made from transparent PVC
- Straight tip
- Available with vent, without vent, or with on/off switch
- Single and double wrap options
- Peel open pouch
- Single patient use


PRODUCT DATA

- Length: 27cm
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Material: PVC (Vent: ABS, PE)
- Shelf Life: 3 Years
- ARTG Number: 224143
- MRI and CT suitable

STANDARDS

- ISO 10993-12:2012 Biological evaluation of medical devices — Part 12: Sample preparation and reference materials
- EN ISO11607-1:2009 Packaging for terminally sterilized medical devices — Part 1: Requirements for materials, sterile barrier systems and packaging systems
- EN 556-1:2001 - Sterilization of medical devices

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Straight Tip
2	No Vent
3	Vent
4	Stepped Connector
5	On/Off Switch

Yankauer Suction Handles Straight Tip

Yankauer Handle - Straight Tip				
RE-ORDER	DESCRIPTION	Type	SIZE	CARTON/QTY
AN170000	Yankauer Handle - Straight Tip without Vent	Single Wrap	8mm	50
AN170000NS	Yankauer Handle - Straight Tip Non-Sterile	Single Wrap	8mm	50
AN170001	Yankauer Handle - Straight Tip with Vent	Single Wrap	8mm	50
AN170005	Yankauer Handle - Straight Tip with On/Off Switch	Single Wrap	8mm	50
AN174000	Yankauer Handle - Straight Tip without Vent	Double Wrap	8mm	50
AN174001	Yankauer Handle - Straight Tip with Vent	Double Wrap	8mm	50
AN174005	Yankauer Handle - Straight Tip with On/Off Switch	Double Wrap	8mm	50