

Yankauer Suction Handles Crown Tip

DESCRIPTION

m|devices Yankauer Suction Handles consists of suction tube. The suction tips are available as a crown tip offering with vent, without vent, or with an on/off switch. Smooth and clean edges for atraumatic suctioning made from a transparent medical grade PVC

FEATURES & BENEFITS

- Made from transparent PVC
- Crown tip
- Available with vent, without vent, or with on/off switch
- Single and double wrap options
- Peel open pouch
- Single patient use


PRODUCT DATA

- Length: 27cm
- Not Made with DEHP
- Not Made with Natural Rubber Latex
- Material: PVC (Vent: ABS, PE)
- Shelf Life: 3 Years
- ARTG Number: 224143
- MRI and CT suitable

STANDARDS

- ISO 10993-1:2009 Biological evaluation of medical devices
- ISO 11737-2:2009(en) Sterilisation of medical devices — Microbiological methods — Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilisation process
- ISO 14971:2019 Application of risk management to medical devices

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Crown Tip
2	No Vent
3	Vent
4	Stepped Connector
5	On/Off Switch

Yankauer Suction Handles Crown Tip

Yankauer Handle - Crown Tip				
RE-ORDER	DESCRIPTION	Type	SIZE	CARTON/QTY
AN170002	Yankauer Handle - Crown Tip without Vent	Single Wrap	8mm	50
AN170002NS	Yankauer Handle - Crown Tip without Vent Non-Sterile	Single Wrap	8mm	50
AN170003	Yankauer Handle - Crown Tip with Vent	Single Wrap	8mm	50
AN170003NS	Yankauer Handle - Crown Tip with Vent Non-Sterile	Single Wrap	8mm	50
AN170004	Yankauer Handle - Crown Tip with On/Off Switch	Single Wrap	8mm	50
AN174002	Yankauer Handle - Crown Tip without Vent	Double Wrap	8mm	50
AN174003	Yankauer Handle - Crown Tip with Vent	Double Wrap	8mm	50
AN174004	Yankauer Handle - Crown Tip with On/Off Switch	Double Wrap	8mm	50