

Yankauer Suction Handles Paediatric

DESCRIPTION

m|devices Paediatric Yankauer Suction Handles consists of suction tip and suction tube. The suction tips are available in a regular or fine tip, with or without vent. Smooth and clean edges for atraumatic suctioning and made from a transparent medical grade PVC.

FEATURES & BENEFITS

- Made from transparent PVC
- Regular or fine tip
- Available with vent, without vent
- Single and double wrap options
- Peel open pouch
- Single patient use


PRODUCT DATA

- Length: 27cm
- Outer diameter: 8mm
- Regular Tip OD: 6mm
- Fine Tip OD: 4mm
- Not Made with DEHP
- Not Made with Natural Rubber Latex
- Material: PVC (Vent: ABS, PE)
- Shelf Life: 3 Years
- ARTG Number: 224143
- MRI and CT suitable
- Sterile EO

STANDARDS

- ISO 10993-1:2009 Biological evaluation of medical devices
- ISO 11737-2:2009(en) Sterilisation of medical devices — Microbiological methods — Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilisation process
- ISO 14971:2019 Application of risk management to medical devices

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Regular or Fine Tip
2	No Vent
3	Vent

Yankauer Suction Handles Paediatric

Yankauer Handle - Paediatric				
RE-ORDER	DESCRIPTION	Type	SIZE	CARTON/QTY
AN170007	Yankauer Handle - Paediatric Fine Tip with Vent	Single Wrap	4mm	50
AN170008	Yankauer Handle - Paediatric Fine Tip without Vent	Single Wrap	4mm	50
AN170009	Yankauer Handle - Paediatric Regular Tip with Vent	Single Wrap	6mm	50
AN170010	Yankauer Handle - Paediatric Regular Tip without Vent	Single Wrap	6mm	50
AN174007	Yankauer Handle - Paediatric Fine Tip with Vent	Single Wrap	4mm	50
AN174008	Yankauer Handle - Paediatric Fine Tip without Vent	Double Wrap	4mm	50
AN174009	Yankauer Handle - Paediatric Regular Tip with Vent	Double Wrap	6mm	50
AN174010	Yankauer Handle - Paediatric Regular Tip without Vent	Double Wrap	6mm	50