

Swabable

Needleless Access Site

PRODUCT DATA SHEET



DESCRIPTION

m|devices Swabable Needleless Access Site has a split septum and is lipid resistant offering high flow rates of up to 550mL/min and can support up to 300 access events for both fluid injection and fluid aspiration. It is luer activated and designed for use with luer lock and luer slip connectors. No additional devices are required for activation.

FEATURES AND BENEFITS

- Swabable needleless valve
- High flow rates due to unimpeded flow
- Easy multiple access with either luer slip or luer lock syringes
- Low priming volume
- Transparent housing for clear visibility
- Simple design with only three components
- Seven-day microbial challenge with 140 activations maintains integrity
- Flushing study to confirm 100% clearance of blood

PRODUCT DATA

- Flow rate: Up to 550mL/min
- Number of accesses: Up to 300 times
- Priming volume: 0.09 mL
- Pressure rating: 37 Psi
- Lipid resistant: Yes
- Shelf Life: 5 years
- Remarks: DEHP / Latex Free
- Material: Silicon inner & poly carbonate outer
- ARTG Number: 272212
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL

RE-ORDER

HR2452-S

DESCRIPTION

Swabable Needleless Access Site

QUANTITY

600 Pieces/Ctn (12boxes, 50pcs/box)

LOT 888888888

yyyy/mm

yyyy/mm

00 kg

STERILE EO

m|devices

Multigate Medical Devices Pty Ltd
27 Llewellyn Avenue Villawood NSW 2163 Australia
T +61(0)2 8718 2800 F +61(0)2 8718 2801
E info@mdevices.com.au W mdevices.com.au

Multigate Co., Ltd.
27 Llewellyn Avenue Villawood NSW 2163 Australia
T +61(0)2 8718 2800 F +61(0)2 8718 2801
E info@mdevices.com.au W mdevices.com.au

CE 0123

RE-ORDER	SIZE	CARTON QTY
HR2452-S		50