

Infusion Set

PRODUCT DATA SHEET



DESCRIPTION

The Infusion Set includes a vented spike attached to a drip chamber with 15µm filter. The 220cm kink resistant standard tubing with male luer lock rotating collar and is made from medical grade PVC. There is also a roller regulator and a latex free injection site.

FEATURES AND BENEFITS

- Vented perforator
- Drip chamber
- Latex free injection site
- Rigid roller regulator for accurate flow
- Male Luer Lock with rotating collar
- Peel open pouch
- Single use only

PRODUCT DATA

- Tubing length: 220cm
- Tubing diameter: 3mm
- Drip chamber: 20 drops /mL
- Filter: 15µm
- Priming volume: 19.0mL
- Sterile EO
- Latex free
- Material: PVC, PE, ABS, PA, MBS
- DEHP: No
- Shelf Life: 3 Years
- ARTG Number: 176900
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV010000	Infusion Set	220cm	100