



Infusion Set

220cm

PRODUCT DATA SHEET



DESCRIPTION

The Infusion Set includes a vented spike attached to a drip chamber with 15µm filter, roller regulator and 1x NeutralSite™.


FEATURES AND BENEFITS

- Vented perforator
- Male luer lock with rotating collar
- NeutralSite™
- Rigid roller regulator for accurate flow
- Peel open pouch
- Single use only

PRODUCT DATA

- Tubing length: 220cm
- Tubing diameter: 3mm
- Drip chamber: 20 drops /mL
- Filter: 15µm
- Priming volume: 18.0mL
- NeutralSite™: Neutral displacement swabable needleless access site
- Sterile EO
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Not made with Natural Rubber Latex
- Materials: LDPE, ABS, PVC, Nylon filter
- NeutralSite™ Materials: PP, ABS and Silicone
- Shelf Life: 3 Years
- ARTG Number: 304860
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV010003	Infusion Set	220cm	50

m|devices®

m|devices
27 Llewellyn Avenue Villawood
NSW 2163 Australia
T + 61 (0) 2 8718 2800 F + 61 (0) 2 8718 2801

info@mdevices.com.au
mdevices.com.au