

Transfusion Set 220cm

– Injection Site

PRODUCT DATA SHEET



DESCRIPTION

The Gravity Transfusion set includes a hydrophobic vented perforator attached to a drip chamber with a 200µm filter. The kink resistant standard tubing and is made from medical grade PVC. There is a male luer lock with rotating collar, roller regulator and a needle injection site which is 30cm from the distal end.

FEATURES AND BENEFITS

- Hydrophobic vented perforator
- Needle injection site
- Rigid roller regulator for accurate flow
- Male luer lock with rotating collar
- Peel open pouch
- Single use only
- Sterile EO
- Latex free

PRODUCT DATA

- Tubing length: 220cm
- Tubing diameter: 3.0mm
- Drip chamber: 20 drops /mL
- Filter: 200µm
- Priming volume: 19.0mL
- Material: PVC, PE, ABS, PA, MBS
- DEHP: No
- Shelf Life: 5 years
- ARTG Number: 228375
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	SIZE	CARTON QTY	RE-ORDER	SIZE	CARTON QTY
IV020000	220cm	100			