



Transfusion Set

PRODUCT DATA SHEET



DESCRIPTION

The Transfusion Set includes a hydrophobic vented perforator attached to a drip chamber with a 200µm filter. The kink resistant standard tubing is available in two lengths and is made from medical grade PVC. There is a male luer lock with rotating collar, roller regulator and a needleless access site, 30cm from the distal end.


FEATURES AND BENEFITS

- Hydrophobic vented perforator
- Needleless access site
- Rigid roller regulator for accurate flow
- Male luer lock with rotating collar
- Peel open pouch
- Single use only

PRODUCT DATA

- Tubing length: 405cm
- Tubing diameter: 3.0mm
- Drip chamber: 20 drops /mL
- Filter: 200µm
- Priming volume: 25.0mL
- Sterile EO
- Latex free
- Material: PVC, PE, ABS, PA, MBS
- DEHP: No
- Shelf Life: 3 Years
- ARTG Number: 228375
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



| RE-ORDER | DESCRIPTION | SIZE | CARTON/QTY |
|----------|-----------------|-------|------------|
| IV020002 | Transfusion Set | 405cm | 50 |