

# Transfusion Set with NRV and 2x Needleless Access Site

## PRODUCT DATA SHEET



### DESCRIPTION

The Gravity Transfusion Set includes a single vented spike attached to a drip chamber with a 200µm filter. The 245cm kink resistant standard tubing with male luer lock connector and rotating collar. It is made from medical grade PVC. There are two swabable needleless access sites and a non-return valve.


### FEATURES AND BENEFITS

- Vented perforator
- Rigid roller regulator for accurate flow
- 2 x swabable needleless access sites
- Non-Return Valve
- Click clamp
- Male luer lock with rotating collar
- Peel open pouch
- Single use only
- Sterile EO
- Latex free

### PRODUCT DATA

- Tubing length: 245cm
- Tubing diameter: 3mm
- Drip chamber: 20 drops /mL
- Filter: 200µm
- Priming volume: 22.0mL
- Material: PVC, PE, ABS, PA, MBS
- DEHP: No
- Shelf Life: 5 years
- ARTG Number: 228375
- The product is MRI and CT suitable

### PRODUCT PACKAGING LABEL



RE-ORDER	SIZE	CARTON QTY
IV020003	245cm	100