

Transfusion Set

with 6cm Side Arm
251cm

DESCRIPTION

The Transfusion Set with Side Arm includes a single spike attached to a drip chamber with a 200µm filter. The 251cm kink resistant standard tubing with male luer lock rotating collar is made from medical grade PVC. There are two NeutralSite™, two non-return valves and a 6cm side arm.

FEATURES & BENEFITS

- Vented perforator
- Rigid roller regulator for accurate flow
- 2 x NeutralSite™
- 6cm Side arm
- Lipid resistant
- 2 x Non-return valves
- Male luer lock with rotating collar
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter 3.0mm
- Tubing length 251cm
- Priming volume: 24.0mL
- Drip chamber: 20 drops /mL
- Filter: 200µm
- NeutralSite™: Neutral displacement swabable needleless access site
- Sterile EO
- ARTG Number: 354953
- The product is MRI and CT suitable
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Not made with Natural Rubber Latex
- Material: PC, PVC, PE, ABS, PP
- NeutralSite™ Materials: PP, ABS and Silicone
- Shelf Life: 3 Years

ISO STANDARDS

- EN ISO 8536-4:2013 Infusion equipment for medical use- part 4: infusion set
- ISO 594/1; Conical Fittings with a 6% (Luer) Taper for Syringes, Needles and certain other Medical Equipment – Part 1: General Requirements
- ISO 594/2; Conical Fittings with a 6% (Luer) Taper for Syringes, Needles and certain other Medical Equipment – Part 2: Lock Fittings
- ISO14971:2016 Medical Devices- Application of risk management to medical device
- ISO13485:2016 Medical Devices- Quality management systems – Requirements for regulatory purposes

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Drip Chamber
2	Back Check valve
3	Clamp
4	NeutralSite™ x 2
5	Roller Clamp
6	Tubing
7	Male Luer Lock

Transfusion Set			STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV020010	Transfusion Set with 6cm Side Arm	251cm	100