



Transfusion Pump Set 220cm

PRODUCT DATA SHEET



DESCRIPTION

The Gravity Transfusion Pump Set includes a single vented spike attached to a flexible drip chamber with a 200µm filter. The 220cm kink resistant tubing with male luer lock and rotating collar is made from medical grade PVC. There is a manual hand pump, roller regulator and swabable needleless access valve.

FEATURES AND BENEFITS

- Vented perforator
- Manual hand pump
- Rigid roller regulator for accurate flow
- Swabable needleless access valve
- Male luer lock with rotating collar
- Peel open pouch
- Single use only
- Sterile EO
- Latex free

PRODUCT DATA

- Tubing length: 220cm
- Tubing diameter: 3mm
- Drip chamber: 20 drops/mL
- Filter: 200µm
- Priming volume: 67.0mL suitable
- Material: PVC, PE, ABS, PA, MBS
- DEHP: No
- Shelf Life: 5 years
- ARTG Number: 228375
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	SIZE	CARTON QTY
IV030001	220cm	50