



Transfusion Set with Manual Pump

PRODUCT DATA SHEET



DESCRIPTION

The Transfusion Set includes a single vented spike attached to a flexible drip chamber with 200µm filter. The 265cm kink resistant tubing, with male luer lock connector and rotating collar, is made from medical grade PVC. There is a manual hand pump, non-return valve, roller clamp and two swabable needleless access sites.


FEATURES AND BENEFITS

- Vented perforator
- Manual hand pump
- Rigid roller regulator for accurate flow
- 2 x Swabable needleless access valve
- Male luer lock with rotating collar
- Non-return valve
- Peel open pouch
- Single use only

PRODUCT DATA

- Tubing length: 265cm
- Tubing diameter: 3.0mm
- Drip chamber: 20 drops /mL
- Filter: 200µm
- Priming volume: 74.0mL
- Sterile EO
- Latex free
- Material: PVC, PE, ABS, PA, MBS
- DEHP: No
- Shelf Life: 5 Years
- ARTG Number: 228375
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QT
IV030004	Transfusion Set with Manual Pump	265cm	50



m|devices
27 Llewellyn Avenue Villawood
NSW 2163 Australia
T + 61 (0) 2 8718 2800 F + 61 (0) 2 8718 2801

info@mdevices.com.au
mdevices.com.au