


Transfusion Pump Set with Luer Break

PRODUCT DATA SHEET



DESCRIPTION

The Transfusion Pump Set with Luer Break includes an integrated Needleless Access Site and manual hand pump. The 300cm kink resistant tubing, with male luer lock with rotating collar connector, is made from medical grade PVC.


FEATURES AND BENEFITS

- Male luer lock with rotating collar
- Swabable Needleless Access Site
- Luer break with female luer
- Non-Return Valve
- Roller regulator for accurate flow x3
- Rapid infusion hand pump
- Non-vented double spike perforators
- Flexible double drip chamber
- Peel open pouch

PRODUCT DATA

- Tubing length: 300cm
- Tubing diameter: 3mm ID
- Drip chamber: 20 drops /mL
- Filter: 200 micron filter
- Sterile EO
- Latex free
- Material: Medical grade PVC, PP, ABS, silicone
- Shelf Life: 5 Years
- ARTG Number: 228375

PRODUCT PACKAGING LABEL



RE-ORDER

DESCRIPTION

SIZE

CARTON/QTY

IV030008

Transfusion Pump Set with Luer Break

300cm

50