

# Burette Set

## 150mL In-Line Port



### PRODUCT DATA SHEET

#### DESCRIPTION

The m|devices Burette Set In-Line Port has a 150mL capacity vented chamber with a floating safety shut off valve. There is also a vented spike, roller regulator and a Swabable Needleless Access Site at one end and In-line Port with click clamp at the other.

#### FEATURES AND BENEFITS

- Vented perforator
- 150mL Burette
- Rigid roller regulator for accurate flow
- Swabable Needleless Access Site
- Floating safety shut off valve
- Click clamp
- In-line Port
- Peel open pouch
- Single use only

#### PRODUCT DATA

- Burette capacity: 150mL
- Tubing length: 25cm
- Drip chamber: 20 drops /mL
- Priming volume: 1.8mL
- Sterile EO
- Latex free
- Material: PVC, PE, ABS
- DEHP: Yes
- Shelf Life: 5 Years
- ARTG Number: 169577
- The product is MRI and CT suitable
- Not to be used with blood or blood products

#### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV040001	Burette Set In-Line Port	25cm	40