

Burette with Needleless Luer Access

PRODUCT DATA SHEET



DESCRIPTION

The Burette with Needleless Luer Access has a 150mL capacity vented drip chamber at 60 drops/mL, with a floating safety shut off valve. There is a vented spike, roller regulator and a swabable needleless access site at one end. While the other has a 195cm tubing with an additional swabable valve and roller regulator, click clamp and a male luer lock connector with rotating collar.

FEATURES AND BENEFITS

- Vented Perforator
- 2 x Rigid Roller Regulator for accurate flow
- 2 x Swabable valves
- 150mL Burette
- Floating safety shut off valve
- Tubing with click clamp
- Male luer lock with rotating collar
- Peel open pouch
- Single use only
- Sterile EO
- Latex free

PRODUCT DATA

- Burette capacity: 150mL
- Tubing length: 195cm
- Tubing diameter: 3mm
- Drip chamber: 60 drops/mL
- Total length: 255cm
- Priming volume: 21.0mL
- Material: PVC, PE, ABS, PA, MBS
- DEHP: Yes
- Shelf Life: 5 Years
- ARTG Number: 169577
- The product is MRI and CT suitable
- Not to be used with blood or blood products

PRODUCT PACKAGING LABEL



RE-ORDER	SIZE	CARTON QTY
----------	------	------------

IV040002	230cm	25
----------	-------	----