

Microbore Extension Set 25cm

DESCRIPTION

The Microbore Extension Set is a microbore tubing with female to male luer lock ends and rotating collar, which is 25cm in length.

FEATURES & BENEFITS

- Female to male luer lock
- Rotating collar
- Microbore tubing
- Transparent, soft, smooth tubing
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter ID 1.0mm OD 2.1mm
- Tubing length 25cm
- Priming volume: 0.20mL
- Sterile EO
- ARTG Number: 228230
- The product is MRI and CT suitable
- Not Made with Natural Rubber Latex
- Not made with DEHP
- Materials: PE, DEHP Free PVC, ABS, PC
- Shelf Life: 3 Years

ISO STANDARDS

- EN ISO 8536-4 Infusion equipment for medical use- Part 4: Infusion sets for single use, gravity feed
- EN ISO 594-1 Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment
—Part 1: General requirements
- EN ISO 594-2 Conical fittings with 6% (Luer) taper for syringes, needles and certain other medical equipment
—Part 2: Lock fittings

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Female Luer Lock
3	Male Luer Lock
4	Cap
5	Tubing

Microbore Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV050000	Microbore Extension Set	25cm	100