

Microbore Extension Set 150cm

DESCRIPTION

The microbore extension set is a microbore tubing with female to male luer lock ends and rotating collar, this set is 150cm in length.

FEATURES & BENEFITS

- Female to male luer lock
- Rotating collar
- Microbore tubing
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter ID 1.0mm OD 2.1mm
- Tubing length 150cm
- Priming volume: 1.2mL
- Sterile EO
- ARTG Number: 354371
- The product is MRI and CT suitable
- Not Made with Natural Rubber Latex
- Not made with DEHP
- Materials: DEHP Free PVC, PE, PC, ABS
- Shelf Life: 3 Years

ISO STANDARDS

- EN ISO 8536-4 Infusion equipment for medical use- Part 4:Infusion sets for single use, gravity feed
- EN ISO 594-1 Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment
—Part 1: General requirements
- EN ISO 594-2 Conical fittings with 6 % (Luer) taper for syringes, needles and certain other medical equipment
—Part 2:Lock fittings

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Female Luer Lock
3	Male Luer Lock
4	Cap
5	Tubing

Microbore Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV050002	Microbore Extension Set	150cm	100