

Microbore Extension Set 200cm

DESCRIPTION

The Microbore Extension Set is a microbore tubing with female to male luer lock ends and rotating collar, which is 200cm in length.

FEATURES & BENEFITS

- Female to male luer lock
- Rotating collar
- Microbore tubing
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter ID 1.0mm OD 2.1mm
- Tubing length 200cm
- Priming volume: 1.6mL
- Sterile EO
- ARTG Number: 354371
- The product is MRI and CT suitable
- Not Made with Natural Rubber Latex
- DEHP NO = Not detected [Result less than the Method Detection Limit EN 14372]
- Materials: DEHP Free PVC, PC, PE, ABS
- Shelf Life: 3 Years

ISO STANDARDS

- EN ISO 8536-4 Infusion equipment for medical use- Part 4:Infusion sets for single use, gravity feed
- EN ISO 594-1 Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment
—Part 1: General requirements
- EN ISO 594-2 Conical fittings with 6 % (Luer) taper for syringes, needles and certain other medical equipment
—Part 2:Lock fittings

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Female Luer Lock
3	Male Luer Lock
4	Cap
5	Tubing

Microbore Extension Set

STERILE

RE-ORDER

DESCRIPTION

LENGTH

CARTON/QTY

IV050003

Microbore Extension Set

200cm

100