


# Microbore Extension Set

200cm

## PRODUCT DATA SHEET



### DESCRIPTION

The Microbore Extension Set is a microbore tubing with female to male luer lock ends and rotating collar, which is 200cm in length.

### FEATURES AND BENEFITS

- Female to male luer lock
- Rotating collar
- Microbore tubing
- Latex free
- Peel open pouch
- Single use only

### PRODUCT DATA

- Tubing diameter 1mm
- Tubing length 200cm
- Priming volume: 1.57mL
- 290psi rated
- Sterile EO
- Material: PVC
- Shelf Life: 3 Years
- ARTG Number: 228230
- The product is MRI and CT suitable
- DEHP free

### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV050003	Microbore Extension Set	200cm	100