

# Microbore Extension Set 15cm

## DESCRIPTION

The m|devices Microbore Extension Set is 15cm in length with a male luer lock and rotating collar to a swabable NeutralSite™ with a removable slide clamp.

## FEATURES & BENEFITS

- Male Luer Lock with rotating collar
- Neutralsite™ with high flow valves
- Removable slide clamp
- Peel open pouch
- Single use only


## PRODUCT DATA

- Tubing inner diameter: 1.0mm
- Tubing length 15cm
- Priming volume: 0.2mL
- Forward pressure resistance 400psi
- NeutralSite™ needleless access site
- Sterile EO
- ARTG Number: 354371
- DEHP NO = Not detected [Result less than the Method Detection Limit EN 14372]
- Not made with Natural Rubber Latex
- Material: ABS, Silicone, PP, DEHP free PVC, PE
- Shelf Life: 3 Years

## ISO STANDARDS

- EN ISO 8536-4 Infusion equipment for medical use- part 4: infusion set
- ISO 80639-7 Small-bore connectors for liquids and gases in healthcare applications – Part 7: Connectors for intravascular or hypodermic applications
- ISO 14971 Medical Devices- Application of risk management to medical device
- ISO 13485 Medical Devices- Quality management systems – Requirements for regulatory purposes

## PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	NeutralSite™
2	Tubing
3	Slide Clamp
4	Male Luer Lock
5	Cap

## Infusion Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV050015	Microbore Extension Set	15cm	200