



Microbore Extension Set

30cm

PRODUCT DATA SHEET



DESCRIPTION

The Microbore Extension Set has a 1mm microbore tubing with female to male luer lock fixed end, slide clamp and is 30cm in length.

FEATURES AND BENEFITS

- Female to fixed male luer lock
- Slide clamp
- Microbore tubing
- Latex free
- Peel open pouch
- Single use only

PRODUCT DATA

- Tubing diameter: 1mm OD: 2.05mm
- Tubing length 30cm
- Priming volume: 0.023mL
- Sterile EO
- DEHP free
- Material: PVC
- Shelf Life: 3 Years
- ARTG Number: 228230
- The product is MRI and CT suitable
- Lipid resistant: Yes

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV050030	Microbore Extension Set	30cm	100

m|devices[®]

m|devices
27 Llewellyn Avenue Villawood
NSW 2163 Australia
T + 61 (0) 2 8718 2800 F + 61 (0) 2 8718 2801

info@mdevices.com.au
mdevices.com.au