

Microbore Extension Set

30cm

DESCRIPTION

The Microbore Extension Set has a 1mm microbore tubing with female to male luer lock fixed end, slide clamp and is 30cm in length.

FEATURES & BENEFITS

- Female to fixed male luer lock
- Slide clamp
- Microbore tubing
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter: 1mm OD: 2.05mm
- Tubing length 30cm
- Priming volume: 0.023mL
- Sterile EO
- The product is MRI and CT suitable
- Lipid resistant
- DEHP NO = Not detected [Result less than the Method Detection Limit EN 14372]
- Not made with Natural Rubber Latex
- Material: PVC
- Shelf Life: 3 Years
- ARTG Number: 354371

ISO STANDARDS

- EN ISO 8536-4 Infusion equipment for medical use- part 4: infusion set
- ISO 80639-7 Small-bore connectors for liquids and gases in healthcare applications – Part 7: Connectors for intravascular or hypodermic applications
- ISO 14971 Medical Devices- Application of risk management to medical device

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Female Luer Lock
3	Tubing
4	Slide Clamp
5	Male Luer Lock
6	Cap

Microbore Extension Set

STERILE

RE-ORDER

DESCRIPTION

LENGTH

CARTON/QTY

IV050030

Microbore Extension Set

30cm

100