



Microbore Extension Set

140cm

PRODUCT DATA SHEET



DESCRIPTION

The Microbore Extension Set has a 1mm microbore tubing with female to male luer lock fixed end, slide clamp and is 140cm in length.

FEATURES AND BENEFITS

- Female to fixed male luer lock
- Slide clamp
- Microbore tubing
- Latex free
- Peel open pouch
- Single use only

PRODUCT DATA

- Tubing diameter 1mm OD: 2.05mm
- Tubing length 140cm
- Priming volume: 1.1mL
- Sterile EO
- DEHP free
- Material: PVC
- Shelf Life: 3 Years
- ARTG Number: 228230
- The product is MRI and CT suitable
- Lipid resistant: Yes

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV050140	Microbore Extension Set	140cm	100