



# Microbore Extension Set

180cm

## PRODUCT DATA SHEET



### DESCRIPTION

The Microbore Extension Set has a 1mm microbore tubing with female to male luer lock fixed end, slide clamp and is 180cm in length.


### FEATURES AND BENEFITS

- Female to fixed male luer lock
- Slide clamp
- Microbore tubing
- Latex free
- Peel open pouch
- Single use only

### PRODUCT DATA

- Tubing diameter: 1mm OD: 2.05mm
- Tubing length 180cm
- Priming volume: 1.5mL
- Sterile EO
- DEHP free
- Material: PVC
- Shelf Life: 3 Years
- ARTG Number: 228230
- The product is MRI and CT suitable
- Lipid resistant: Yes

### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV050180	Microbore Extension Set	180cm	100