

Microbore Extension Set

140cm

DESCRIPTION

The Microbore Extension Set has a 1mm microbore tubing with female to male luer lock fixed end, slide clamp and is 180cm in length.

FEATURES & BENEFITS

- Female to fixed male luer lock
- Slide clamp
- Microbore tubing
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter: ID-1.0mm, O.D-2mm
- Tubing length 180cm
- Priming volume: 1.5mL
- Sterile EO
- The product is MRI and CT suitable
- Lipid resistant
- Not made with DEHP
- Not made with Natural Rubber Latex
- Material: PVC
- Shelf Life: 3 Years
- ARTG Number: 354371

ISO STANDARDS

- ISO 13485-2016 Quality management systems. Requirements for regulatory purposes
- ISO 1497-2012 Medical devices. Application of risk management to medical devices
- ISO 11135-2014 Sterilisation of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilisation process for medical devices

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Female Luer Lock
3	Tubing
4	Slide Clamp
5	Male Luer Lock
6	Cap

Microbore Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV050180	Microbore Extension Set	180cm	100