




Minimum Volume Extension Set

600cm

PRODUCT DATA SHEET



DESCRIPTION

The Minimum Volume Extension Set is a 6m tubing female luer lock to male luer lock ends and rotating collar. The long 6m length is designed for MRI use.


FEATURES AND BENEFITS

- Female to fixed male luer lock with rotating collar
- Long 6m tubing for MRI
- Latex free
- Peel open pouch

PRODUCT DATA

- Tubing diameter: 1.5mm
- Tubing length: 600cm
- Priming volume: 10.6mL
- Sterile EO
- DEHP free
- Material: PVC
- Shelf Life: 3 Years
- ARTG Number: 228230
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV051002	Minimum Volume Extension Set	600cm	50

m|devices

m|devices
27 Llewellyn Avenue Villawood
NSW 2163 Australia
T + 61 (0) 2 8718 2800 F + 61 (0) 2 8718 2801

info@mdevices.com.au
mdevices.com.au