

## Minimum Volume Extension Set

### 600cm

#### DESCRIPTION

The Minimum Volume Extension Set is a 6m tubing female luer lock to male luer lock ends and rotating collar. The long 6m length is designed for MRI use.

#### FEATURES & BENEFITS

- Female to fixed male luer lock with rotating collar
- Long 6m tubing for MRI
- Peel open pouch


#### PRODUCT DATA

- Tubing diameter: I.D- 1.5mm, O.D- 3mm
- Tubing length: 600cm
- Priming volume: 10.6mL
- Sterile EO
- The product is MRI and CT suitable
- ARTG Number: 354371
- Not made with DEHP
- Not made with Natural Rubber Latex
- Material: PVC
- Shelf Life: 3 Years

#### ISO STANDARDS

- ISO 13485-2016 Quality management systems. Requirements for regulatory purposes
- ISO 1497-2012 Medical devices. Application of risk management to medical devices
- ISO 11135-2014 Sterilisation of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilisation process for medical devices

#### PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Female Luer Lock
3	Male Luer Loc with Rotating Collar
4	Cap - Vented, FLL suit MLL
5	Tube

#### Minimum Volume Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV051002	Minimum Volume Extension Set	600cm	50