

# Minimum Volume Extension Set

## 165cm

### DESCRIPTION

The Minimum Volume Extension Set is a 165cm tubing with female luer lock to male luer lock ends with rotating collar. It includes a click clamp and a needle injection site.

### FEATURES & BENEFITS

- Female to fixed male luer lock with rotating collar
- Click clamp
- Needle injection site
- Microbore tubing
- Peel open pouch
- Single use only


### PRODUCT DATA

- Tubing diameter: I.D- 1.5mm, O.D- 3mm
- Tubing length: 165cm
- Priming volume: 3mL
- Sterile EO
- The product is MRI and CT suitable
- ARTG Number: 354371
- Not made with DEHP
- Not made with Natural Rubber Latex
- Material: PVC
- Shelf Life: 3 Years

### ISO STANDARDS

- ISO 13485-2016 Quality management systems. Requirements for regulatory purposes
- ISO 1497-2012 Medical devices. Application of risk management to medical devices
- ISO 11135-2014 Sterilisation of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilisation process for medical devices

### PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Male Luer Lock
3	Tubing
4	Y-Injection Site
5	Clamp
6	Female Luer Lock
7	Cap

### Minimum Volume Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV051003	Minimum Volume Extension Set	165cm	100