



Minimum Volume Extension Set

165cm

PRODUCT DATA SHEET



DESCRIPTION

The Minimum Volume Extension Set is a 165cm tubing with female luer lock to male luer lock ends with rotating collar. It includes a click clamp and a needle injection site.



FEATURES AND BENEFITS

- Female to fixed male luer lock with rotating collar
- Click clamp
- Needle injection site
- Microbore tubing
- Latex free
- Peel open pouch
- Single use only

PRODUCT DATA

- Tubing diameter: 1.5mm
- Tubing length: 165cm
- Priming volume: 3mL
- Sterile EO
- DEHP free
- Material: PVC
- Shelf Life: 3 Years
- ARTG Number: 228230
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV051003	Minimum Volume Extension Set	165cm	100



m|devices
 27 Llewellyn Avenue Villawood
 NSW 2163 Australia
 T + 61 (0)2 8718 2800 F + 61 (0)2 8718 2801

info@mdevices.com.au
 mdevices.com.au