




Minimum Volume Extension Set

165cm

PRODUCT DATA SHEET



DESCRIPTION

The Minimum Volume Extension Set is a microbore tubing with female to male luer lock ends and rotating collar. There is a variety of lengths available. The Minimum Volume Extension Set is a 165cm tubing with female to male luer lock ends with rotating collar. It includes a click clamp and a swabable Needleless Access Site.

FEATURES AND BENEFITS

- Female to fixed male luer lock with rotating collar
- Click clamp
- Needleless Access Site
- Latex free
- Peel open pouch
- Single use only

PRODUCT DATA

- Tubing diameter: 1.5mm
- Tubing length: 165cm
- Priming volume: 3mL
- Sterile EO
- DEHP free
- Material: PVC
- Shelf Life: 3 Years
- ARTG Number: 228230
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QT
IV051004	Minimum Volume Extension Set	165cm	100

m|devices®

m|devices
27 Llewellyn Avenue Villawood
NSW 2163 Australia
T + 61 (0) 2 8718 2800 F + 61 (0) 2 8718 2801

info@mdevices.com.au
mdevices.com.au