



Minimum Volume Extension Set

20cm

PRODUCT DATA SHEET



DESCRIPTION

The Minimum Volume Extension Set is a 20cm tubing with female luer lock to male luer lock and rotating collar. It includes a back-check valve and a swabable Needleless Access Site.


FEATURES AND BENEFITS

- Female to male luer lock with rotating collar
- Microbore tubing
- Back-check valve
- Needleless Access Site
- Single use only
- Peel open pouch

PRODUCT DATA

- Tubing diameter: 1.5mm
- Tubing length: 20cm
- Priming volume: 0.35mL
- Sterile EO
- DEHP: No
- Latex free
- Material: PVC
- Shelf Life: 3 Years
- ARTG Number: 228230
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV051006	Minimum Volume Extension Set	20cm	100



m|devices
 27 Llewellyn Avenue Villawood
 NSW 2163 Australia
 T + 61 (0) 2 8718 2800 F + 61 (0) 2 8718 2801

info@mdevices.com.au
 mdevices.com.au