




Standard Bore Extension Set

75cm

PRODUCT DATA SHEET



DESCRIPTION

The Standard Bore Extension Set (75cm) is a female to male luer lock line with rotating collar and a slide clamp.


FEATURES AND BENEFITS

- Female to male luer lock with rotating collar
- Slide clamp
- Peel open pouch
- Single use only
- Sterile EO
- DEHP free
- Latex free

PRODUCT DATA

- Tubing length: 75cm
- Tubing diameter: 3mm
- Priming volume (75cm): 5.31mL
- Priming volume (150cm): 10.6mL
- Material: PVC, PE, PC
- Shelf Life: 5 Years
- ARTG Number: 228230
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	SIZE	CARTON QTY
IV052000	75cm	100