

Standard Bore Extension Set 75cm

DESCRIPTION

The Standard Bore Extension Set (75cm) is a female to male luer lock line with rotating collar and a slide clamp.

FEATURES & BENEFITS

- Female to male luer lock with rotating collar
- Slide clamp
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter: I.D- 3mm O.D- 4mm
- Tubing length: 75cm
- Priming volume (75cm): 5.31mL
- Priming volume (150cm): 10.6mL
- Sterile EO
- The product is MRI and CT suitable
- Lipid resistant: Yes
- Not made with DEHP
- Not made with Natural Rubber Latex
- Material: PVC, PE, PC
- Shelf Life: 3 Years
- ARTG Number: 354371

ISO STANDARDS

- ISO 13485-2016 Quality management systems. Requirements for regulatory purposes
- ISO 1497-2012 Medical devices. Application of risk management to medical devices
- ISO 11135-2014 Sterilisation of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilisation process for medical devices

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Female Luer Lock
2	Slide Clamp
3	Tubing
4	Male Luer Lock with RC
5	Cap

Standard Bore Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV052000	Standard Bore Extension Set	75cm	100