




Standard Bore Extension Set

150cm

PRODUCT DATA SHEET



DESCRIPTION

The Standard Bore Extension Set 150cm has a female to male luer lock with rotating collar distal connector and a slide clamp.


FEATURES AND BENEFITS

- Female to male luer lock with rotating collar
- Slide clamp
- Latex free
- Peel open pouch
- Single use only

PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 150cm
- Priming volume (150cm): 10.6mL
- DEHP free
- Sterile EO
- Material: PVC, PE, PC
- Shelf Life: 3 Years
- ARTG Number: 228230
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV052001	Standard Bore Extension Set	150cm	100



m|devices
 27 Llewellyn Avenue Villawood
 NSW 2163 Australia
 T + 61 (0)2 8718 2800 F + 61 (0)2 8718 2801

info@mdevices.com.au
 mdevices.com.au