

Standard Extension Set

495cm



PRODUCT DATA SHEET

DESCRIPTION

The m|devices Standard Extension Set is 495cm in length with a standard injection site and a female to male luer lock line with rotating collar (RC).


FEATURES AND BENEFITS

- Standard injection site
- Standard bore tubing
- Female luer lock to male lock end with RC
- Made from non-toxic medical grade PVC
- Transparent, soft, smooth tubing
- Peel open pouch
- Single use only
- Latex free

PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 495cm
- Priming volume: 34.99mL
- DEHP free
- Sterile EO
- Material: PVC, PE, PC
- Shelf Life: 3 Years
- ARTG Number: 228230
- The product is MRI and CT suitable
- Lipid resistant: No

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV052100	Standard Extension Set	495cm	50